


Product Datasheet

ATP5L antibody (orb517860)

Description	ATP5L antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human ATP synthase subunit g, mitochondrial protein (2-103AA)
Target	ATP5L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075964
Dilution Range	IHC-P:1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...