


## Product Datasheet

### Collagen IV antibody (orb517876)

|                            |                                                                    |
|----------------------------|--------------------------------------------------------------------|
| <b>Description</b>         | Collagen IV antibody                                               |
| <b>Species/Host</b>        | Rabbit                                                             |
| <b>Reactivity</b>          | Human                                                              |
| <b>Conjugation</b>         | Unconjugated                                                       |
| <b>Tested Applications</b> | ELISA, IHC                                                         |
| <b>Immunogen</b>           | Recombinant Human Collagen alpha-5(IV) chain protein (1184-1235AA) |
| <b>Target</b>              | COL4A5                                                             |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol        |
| <b>Form/Appearance</b>     | Liquid                                                             |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.      |
| <b>Note</b>                | For research use only                                              |
| <b>Isotype</b>             | IgG                                                                |
| <b>Clonality</b>           | Polyclonal                                                         |
| <b>Uniprot ID</b>          | <b>P29400</b>                                                      |
| <b>Dilution Range</b>      | IHC-P:1:20-1:200                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                    |



Immunohistochemistry  
of paraffin-embedde...