

Product Datasheet

HOXC6 antibody (HRP) (orb517897)

Description	HOXC6 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Homeobox protein Hox-C6 protein (1-120AA)
Target	HOXC6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09630
Expiration Date	12 months from date of receipt.