


Product Datasheet

NRL antibody (orb517932)

Description	NRL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Neural retina-specific leucine zipper protein (24-124AA)
Target	NRL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54845
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...