

Product Datasheet

ERK1/2 (phospho-Tyr204) antibody (orb5180)

Description	ERK1/2 (phospho-Tyr204) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat MAPK3 around the phosphorylation site of Tyr204 (TE(p-Y)VA)
Target	phospho-MAPK3 (Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q63844
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.