


Product Datasheet

CPEB3 antibody (orb518008)

Description	CPEB3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Cytoplasmic polyadenylation element-binding protein 3 protein (414-525AA)
Target	CPEB3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NE35
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...