

Product Datasheet

NSL1 antibody (HRP) (orb518017)

Description	NSL1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Kinetochores-associated protein NSL1 homolog protein (225-272AA)
Target	NSL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96IY1
Expiration Date	12 months from date of receipt.