


Product Datasheet

BARD1 antibody (orb518028)

Description	BARD1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human BRCA1-associated RING domain protein 1 protein (428-777AA)
Target	BARD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99728
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...