



Product Datasheet

LSM11 antibody (orb518056)

Description	LSM11 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human U7 snRNA-associated Sm-like protein LSM11 protein (1-141AA)
Target	LSM11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P83369
Dilution Range	WB:1:500-1:5000, IF:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of U251 cell...



Western blot analysis of HepG2 whole ce...