


Product Datasheet

TTC39B antibody (orb518072)

Description	TTC39B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA)
Target	TTC39B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5VTQ0
Expiration Date	12 months from date of receipt.



Western blot
analysis of
HeLa whole
cel...