

Product Datasheet

DPH5 antibody (FITC) (orb518082)

Description	DPH5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Diphthine methyl ester synthase protein (138-235AA)
Target	DPH5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H2P9
Expiration Date	12 months from date of receipt.