



Product Datasheet

LILRB2 antibody (orb518100)

Description	LILRB2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily B member 2 protein (406-563AA)
Target	LILRB2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N423
Dilution Range	WB:1:500-1:5000, IF:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...



Western blot aanalysis of 293T whole cel...