


Product Datasheet

cas3 antibody (orb518108)

Description	cas3 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Thermus thermophilus CRISPR-associated endonuclease/helicase Cas3 protein (120-226AA)
Target	cas3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53VY2
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Recombinant
pr...