

Product Datasheet

SDR antibody (HRP) (orb518121)

Description	SDR antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Short-chain alcohol dehydrogenase protein (1-137AA)
Target	SDR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0A0A0
Expiration Date	12 months from date of receipt.