

Product Datasheet

src antibody (Biotin) (orb518127)

Description	src antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Proto-oncogene tyrosine-protein kinase Src protein (2-146AA)
Target	src
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1JPZ3
Expiration Date	12 months from date of receipt.