

Product Datasheet

MPN_083 antibody (FITC) (orb518138)

Description	MPN_083 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Mycoplasma pneumoniae Uncharacterized lipoprotein MPN_083 protein (22-242AA)
Target	MPN_083
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P75610
Expiration Date	12 months from date of receipt.