

Product Datasheet

CD45 antibody (FITC) (orb518170)

Description	CD45 antibody (FITC)
Species/Host	Rabbit
Reactivity	Rabies virus
Conjugation	FITC
Immunogen	Recombinant Rabies virus Putative uncharacterized protein C protein (1-69AA)
Target	C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15201
Expiration Date	12 months from date of receipt.