


Product Datasheet

zot antibody (orb518180)

Description	zot antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Vibrio cholerae serotype O1 Zona occludens toxin protein (1-399AA)
Target	zot
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38442
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Positive WB
de...