

Product Datasheet

VP3 antibody (FITC) (orb518194)

Description	VP3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Chicken anemia virus Apoptin protein (1-120AA)
Target	VP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54095
Expiration Date	12 months from date of receipt.