


Product Datasheet

vpx antibody (orb518200)

Description	vpx antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx protein (1-112AA)
Target	vpx
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06939
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Recombinant
pr...