

Product Datasheet

pmpD antibody (FITC) (orb518210)

Description	pmpD antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Chlamydia trachomatis Probable outer membrane protein PmpD protein (1244-1531AA)
Target	pmpD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O84818
Expiration Date	12 months from date of receipt.