

Product Datasheet

vcp antibody (HRP) (orb518229)

Description	vcp antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Transitional endoplasmic reticulum ATPase protein (704-806AA)
Target	vcp
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7ZU99
Expiration Date	12 months from date of receipt.