

Product Datasheet

L antibody (FITC) (orb518234)

Description	L antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Vesicular stomatitis Indiana virus RNA-directed RNA polymerase L protein (598-784AA)
Target	L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8B0H0
Expiration Date	12 months from date of receipt.