

Product Datasheet

PG_1424 antibody (FITC) (orb518246)

Description	PG_1424 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Porphyromonas gingivalis Peptidylarginine deiminase protein (44- 556AA)
Target	PG_1424
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9RQJ2
Expiration Date	12 months from date of receipt.