



Product Datasheet

Tubulin beta 3 antibody (orb518249)


Description	Tubulin beta 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fusion protein of Human TUBB3
Target	TUBB3
Preservatives	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13509
Dilution Range	ELISA:1:500-1:5000, WB:1:500-1:1000, IHC-P:1:25-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Human brain gl...



Immunohistochemistry of paraffin-embedde...