

Product Datasheet

CHD5 antibody (FITC) (orb51920)

Description	CHD5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Chromodomain-helicase-DNA-binding protein 5 protein (1600-1699AA)
Target	CHD5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TDI0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.