

Product Datasheet

VPS33A antibody (HRP) (orb51935)

Description	VPS33A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 33A protein (1-300AA)
Target	VPS33A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96AX1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.