


Product Datasheet

Cobra venom factor antibody (orb51978)

Description	Cobra venom factor antibody
Species/Host	Rabbit
Reactivity	Reptile
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Naja kaouthia Cobra venom factor protein (733-984AA)
Target	Cobra venom factor
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q91132
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...