



Product Datasheet

GNB4 antibody (orb52046)


Description	GNB4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Guanine nucleotide-binding protein subunit beta-4 protein (105-185AA)
Target	GNB4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HAV0
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:20-1:200, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat liver(lane ...



Immunohistochemical staining of human pa...



Immunofluorescent analysis of HeLa cells...