

Product Datasheet

ASH1L antibody (FITC) (orb52068)

Description	ASH1L antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Histone-lysine N-methyltransferase ASH1L protein (1906-2160AA)
Target	ASH1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NR48
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.