

Product Datasheet

X antibody (Biotin) (orb52073)

Description	X antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis B virus genotype D subtype ayw Protein X protein (1-154AA)
Target	X
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9QMI3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.