



## Product Datasheet

### GAB2 antibody (orb520882)


|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | GAB2 antibody                                                 |
| <b>Species/Host</b>        | Rabbit                                                        |
| <b>Reactivity</b>          | Human, Mouse, Rat                                             |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                |
| <b>Immunogen</b>           | Fusion protein of Human GAB2                                  |
| <b>Target</b>              | GAB2                                                          |
| <b>Preservatives</b>       | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol       |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>Uniprot ID</b>          | <a href="#">Q9UQC2</a>                                        |
| <b>Dilution Range</b>      | ELISA:1:2000-1:5000, WB:1:500-1:2000, IHC-P:1:50-1:200        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Lane 1-3: A549...



Immunohistochemistry of paraffin-embedde...