

Product Datasheet

MAP1LC3C antibody (Biotin) (orb52185)

Description	MAP1LC3C antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Microtubule-associated proteins 1A/1B light chain 3C protein (1-30AA)
Target	MAP1LC3C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXW4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.