

Product Datasheet

ABCA2 antibody (Biotin) (orb52189)

Description	ABCA2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human ATP-binding cassette sub-family A member 2 protein (1223-1429AA)
Target	ABCA2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZC7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.