

Product Datasheet

Tropomyosin antibody (HRP) (orb52199)

Description	Tropomyosin antibody (HRP)
Species/Host	Rabbit
Reactivity	Crustacean
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Charybdis feriata Tropomyosin protein (1-264AA)
Target	Tropomyosin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9N2R3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.