



Product Datasheet

BAIAP2 antibody (orb52238)


Description	BAIAP2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Brain-specific angiogenesis inhibitor 1-associated protein 2 protein (123-368AA)
Target	BAIAP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UQB8
Dilution Range	WB:1:200-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of Hela whole cell lysate.



Immunofluorescence analysis of HepG2 cells.