

Product Datasheet

BACE2 antibody (HRP) (orb52275)

Description	BACE2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Beta-secretase 2 protein (84-157AA)
Target	BACE2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5Z0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.