

Product Datasheet

APOC2 antibody (orb5228)

Description	APOC2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, WB
Immunogen	KLH conjugated synthetic peptide derived from human APOC2 (23-97/101 aa)
Target	APOC2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9 kDa
Uniprot ID	P02655
Expiration Date	12 months from date of receipt.