


Product Datasheet

ASB4 antibody (orb52306)

Description	ASB4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Ankyrin repeat and SOCS box protein 4 protein (166-321AA)
Target	ASB4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y574
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of ASB4(lane 1), H...