


## Product Datasheet

### ANGPTL2 antibody (orb523068)


|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | ANGPTL2 antibody                                              |
| <b>Species/Host</b>        | Rabbit                                                        |
| <b>Reactivity</b>          | Human, Mouse                                                  |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                |
| <b>Immunogen</b>           | Synthetic peptide of Human ANGPTL2                            |
| <b>Target</b>              | ANGPTL2                                                       |
| <b>Preservatives</b>       | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol       |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>Uniprot ID</b>          | <b>Q9UKU9</b>                                                 |
| <b>Dilution Range</b>      | ELISA:1:1000-1:2000, WB:1:200-1:1000, IHC-P:1:25-1:100        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Mouse intestin...



Immunohistochemistry of paraffin-embedde...