



Product Datasheet

LYVE1 antibody (orb52318)


Description	LYVE1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Lymphatic vessel endothelial hyaluronic acid receptor 1 protein (20-238AA)
Target	LYVE1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5Y7
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver(lan...



Immunohistochemical staining of human to...



Immunofluorescence analysis of MCF-7 cel...