

Product Datasheet

AVPR2 antibody (HRP) (orb52387)

Description	AVPR2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Pig Vasopressin V2 receptor protein (292-370AA)
Target	AVPR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32307
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.