


Product Datasheet

BET1 antibody (orb52398)

Description	BET1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human BET1 homolog protein (1-83AA)
Target	BET1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15155
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...