

Product Datasheet

CEACAM6 antibody (Biotin) (orb52462)

Description	CEACAM6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 6 protein (35-320AA)
Target	CEACAM6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40199
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.