



Product Datasheet

CELF2 antibody (orb52503)


Description	CELF2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human CUGBP Elav-like family member 2 protein (291-381AA)
Target	CELF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95319
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(lane 1), A...



Immunohistochemical staining of human to...



Immunofluorescence analysis of A549 cell...