



## Product Datasheet

### TRIM45 antibody (orb525213)


|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | TRIM45 antibody                                               |
| <b>Species/Host</b>        | Rabbit                                                        |
| <b>Reactivity</b>          | Human, Mouse                                                  |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                |
| <b>Immunogen</b>           | Synthetic peptide of Human TRIM45                             |
| <b>Target</b>              | TRIM45                                                        |
| <b>Preservatives</b>       | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol       |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q9H8W5">Q9H8W5</a>     |
| <b>Dilution Range</b>      | ELISA:1:2000-1:5000, WB:1:500-1:2000, IHC-P:1:25-1:100        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Lane 1-5: 231 ...



Immunohistochemistry of paraffin-embedde...