

Product Datasheet

F12 antibody (Biotin) (orb52534)

Description	F12 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Pig Coagulation factor XII protein (20-371AA)
Target	F12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O97507
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.