

Product Datasheet

GAS7 antibody (FITC) (orb52577)

Description	GAS7 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Growth arrest-specific protein 7 protein (1-153AA)
Target	GAS7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60861
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.