

Product Datasheet

Acetyl Histone H4 antibody (FITC) (orb5258)

Description

Rabbit polyclonal antibody against Acetyl

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLKLH conjugated Synthesised acetylpeptide derived from human Histone H4 around the acetylation site of K16
Target	H4C1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P62805
Dilution Range	Flow-Cyt=1ug/TEST,IF=1:50-200
Expiration Date	12 months from date of receipt.