

Product Datasheet

GNB1 antibody (Biotin) (orb52618)

Description	GNB1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-1 protein (157-248AA)
Target	GNB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62873
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.